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What is the indication for CaReS®-1S?

CaReS®-1S is indicated for a wide range of cartilage defects:

- Partial joint cartilage defects with a damaged cartilage surface without involvement of the subchondral bone plate
- Deep cartilage defects or subchondral joint cartilage defects involving the subchondral bone plate
- Osteochondral defects as Osteochondrosis Dissecans

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Rehabilitation of your patients

Femur condyles and talus role

- Mobilisation on 2nd post-operative day after removal of drainage
- 20 kg partial weight bearing of the operated extremity for 6 weeks, free mobility
- At least 2 x daily for 1 h continuous passive motion
- From the 7th week on load progression 30 kg/2 weeks
- After achievement of full weight bearing cycling and swimming allowed, careful building of the muscles by isometric training
- Jump-, run- and contact sports allowed after 1 year

Patella and glide bearing

- Mobilisation on 2nd post-operative day after removal of drainage. The flexion of the knee is limited by an IROM-splint for 6 weeks
- 20 kg partial weight bearing of the operated extremity for 1 weeks
- From the 2nd week full weight bearing allowed
- From the 6th week increase of the flexion 30°/2 weeks, from 90° flexion the IROM-splint could be removed
- From 90° flexion cycling and swimming allowed, careful building of the muscles by isometric training
- Sports with high risk to fall are allowed after 1 year

CaReS®-1S

Cartilage Regeneration System - One Step

Produktinformation / Product Information

ARTHRO KINETICS
REGENERATING JOINT MOBILITY

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What is CaReS®-1S?

CaReS®-1S is a round implant to treat focal cartilage defects up to 6 cm². CaReS®-1S is a sterile ready-to-use primary cellfree implant based on a dense patented matrix of Collagen type I. CaReS®-1S fills the defect completely, is colonized by cells of the surrounding tissue and leads to the formation of hyaline-like cartilage.

CaReS®-1S is a technology of cartilage regeneration which taps the advantages of an Autologous Chondrocyte Implantation without the need of more than one intervention to get of a biopsy.

Properties of CaReS®-1S:

- Effective reduction of pain and improvement of mobility; improvement of quality of life of YOUR patients
- Proven immigration of cells into the CaReS®-1S matrix and formation of hyaline-like cartilage
- One-time gentle intervention
- Easy to apply and fix; simple and fast OR procedure
- Ready-to-use; CaReS®-1S is at anytime available in your OR; you avoid logistical effort
- No licensing procedure
- CE marked

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

Regenerative Potential von CaReS®-1S:

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
CaReS®-1S - First results:

Results of experimental studies in nude mice and minipigs as well as first clinical results.

Regenerative potential of CaReS®-1S:

- CaReS®-1S remains in the defect and is stable
- CaReS®-1S is colonized by precursor cells migrating from the surrounding tissue (6-12 weeks)
- These cells differentiate into chondrocytes, express Collagen Type II and transform CaReS®-1S into a stable hyaline-like tissue (52 weeks)
- CaReS®-1S integrates very well into the defect and the surrounding tissue. Its regenerative potential can be compared to that of cell-colonized transplants

First results of a clinical study on 18 patients with focal femoral condyle and retropatellar joint defects indicate:

- Evident pain relief
- Successful defect filling
- Formation of hyaline-like cartilage
- Very good compatibility
- Clear tendency towards functional improvement